



ELECTRICAL POWER PLUG – NEMA TYPE

Yamato and other scientific equipment manufacturers do not include electrical power plugs with some equipment, instead it delivers with round terminals at the end of the power cord. The main reason for this practice is the diversity of type of plugs and voltage at the laboratories or facilities around the globe. To help customers to select the proper electrical power plug type, we provide instructions in the equipment user manual.

Particularly for US customers, we recommend NEMA (National Electrical Manufacturers Association) power plugs, most commonly used in the USA, Canada and Japan. The table below shows NEMA electrical power plug type most commonly use on Yamato equipment. Please check the user manual for voltage and amperage. For other countries, please consult NEMA to find the correspondent electrical power plug model.



These models do not include the power plug. Correctly connect the ground to fit with the power supply facility to be connected.

